Designs and performance of three new microprocessor-controlled knee joints

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C-Leg 4, Rheo Knee 3 and Plié 3

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With C-Leg 4 compared to Rheo Knee 3 and Plié 3:

→ **C-Leg 4 compared to Rheo Knee 3**

Increased mean self-selected walking velocity by **0.15 m/s**

Maximum knee flexion in stance phase **3.61° higher**

→ **C-Leg 4 compared to Plié 3**

Microprocessor-controlled adaption of flexion and extension with C-Leg 4 during various gait speeds

→ Closer to physiological gait pattern

Clear superiority due to swing phase control

→ **No automatic flexion or extension damping with Plié 3**

- Characterized by a high correlation between max. knee flexion angle and walking velocity ($R^2 = 0.9$).

→ **All MPKs**

Reliable detection of stance and swing phase

Reliable generation of joint resistance

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

Figure 1: The graph represents the maximum knee angle [°] over walking velocity [m/s]. It is shown that the slope of the C-Leg 4 trend line is closest to the contralateral side, thus reflecting the most natural behaviour.

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Subjects: Four male, unilateral, transfemoral amputees
 Previous prosthesis: Genium
 Amputation causes: not stated
 Mean age: 45.8 yrs (± 12.09 yrs)
 Mean time since amputation: 27.8 yrs (± 10.5 yrs)
 MFCL: K3

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"The C-Leg 4, Plié 3 and Rheo Knee 3 offer a reliable detection of stance and swing phase and the generation of joint resistance to avoid uncontrolled flexion of the knee joint. As shown in the technical analysis, only C-Leg 4 and Rheo Knee 3 allow for microprocessor-controlled adaptation of resistance to different movements and situational requirements. The Plié has manually adjustable elements that generate the resistances. This directly affects the outcome of the biomechanical analysis which examined level walking at different walking velocities. C-Leg 4 showed the most natural knee function when compared to the contralateral side and this was followed by Rheo Knee 3." (Thiele et al., 2018)

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© 2018, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.