## Benefits of microprocessor-controlled prosthetic knees to limited community ambulators: Systematic review

*Journal of Rehabilitation Research & Development (JRRD) 2014; 51(10): 1469-1496.*

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**C-Leg / Compact vs NMPKs** 機械式膝継手

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With C-Leg / Compact compared to NMPKs:

→ **Limited community ambulators (MFCL2) are safer with C-Leg**

Decreased number of stumbles and falls by 80%.

Increased balance in Activities of Daily Living (ADLs)

→ **Limited community ambulators (MFCL2) benefit from C-Leg in walking**

14-25 % faster walking on level ground

20% faster on uneven ground

Improved walking quality on stair and ramp descent

→ **Limited community ambulators (MFCL2) prefer C-Leg over NMPKs**

Up to 90% prefer C-Leg over their previous NMPK

→ **No significant benefits nor trends towards an advantage of NMPKs were reported compared to the C-Leg or Compact**

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
## Outcome measures showing the benefits of C-Leg

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160 outcome measures were analysed in this review. The graph shows the percentage of outcome measures reporting a significant improvement ( $p<0.05$ ) or statistical trend ( $0.05<p<0.1$ ) towards a benefit for C-Leg or if no statistical difference was found. No benefits were reported for the NMPKs when compared to C-Leg.

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
**Subjects:** 57 limited community ambulators  
**Amputation causes:** Dysvascular, PVD or Diabetes (20), Trauma (31), Other (6)  
**Mean age:** 57.1 – 67.1 years  
**MFCL:** MFCL 2

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"The results of this systematic review of clinical trials on interventions with MPKs in individuals with a unilateral TFA and MFCL-2 mobility grade suggest that these subjects may significantly reduce the number of falls and their risk of falling, improve their balance, and better perform activities of community ambulation that are actually categorized as part of the MFCL-3 mobility grade. Because these results have been derived from studies with low to moderate methodological quality in a yet limited number of patients, trial fittings with different types of MPKs (MP stance only or MP stance and swing control) may be considered to evaluate whether an individual benefits from using an MPK compared with NMPKs usually prescribed for MFCL-2 individuals. Criteria for appraising success or failure of the trial fitting based on the 2MWT, AMP, TUG, and ABC have been suggested. Given the challenges to objectify the current general and ambiguous definitions of the MFCLs, an evidence-based and unambiguously quantifiable functional classification or one or more validated outcome measures to corroborate the classification would help better define patient groups to be subjected to clinical research and sharpen coverage and reimbursement criteria." (Kannenberg et al, 2014)

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