Effects of Mobility Grade, Age and Etiology on functional benefit and safety of subjects evaluated in over 1200 C-Leg trial fittings in Germany

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C-Leg and C-Leg Compact

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With C-Leg and Compact compared to the previous prosthesis:

→ **Utilization of functional benefit is high even after one-day test fittings.**

(ranging from 83% to 95%, reduction of walking aids 46%)

→ **Utilization of functional benefit is independent from age, mobility grade and/or etiology.**

Largest correlation < 0.251

→ **Age, mobility grade and/or etiology have no meaningful diagnostic predictive power**

Ranging from 0.7% to 9%

→ **Even within test fitting significant changes in mobility grade can be observed**

50% of subjects increased their mobility grade from MFCL 2 to MFCL 3

22% of subjects increased their mobility grade from MFCL 3 to MFCL 4

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→ Gait was improved by C-Leg

Harmonized gait pattern
(88-95% of subjects)

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
Fear of falling with C-Leg compared to previous prosthesis

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
Subjects: 1223 amputees (Transfemoral 84%, Knee disarticulation 13%, Hip disarticulation 3%)
Previous prosthesis: not reported
Amputation causes: Trauma (43%), vascular disease incl. diabetes (26%), tumor (13%), infection (6%), combat trauma (4%), malformation (2%), others (6%), not specified (13%)
Mean age: 55.6 ± 15.1 yrs
Mean time since amputation: 17.5 ± 18.9 yrs
MFCL: K2 (38.4%), K3 (39.2%), K4 (6.5%) (K1 (6.1%) excluded from analysis)

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Retrospective, cross sectional analysis:

Data of C-Leg or C-Leg Compact trial fittings were gathered from ca. 445 prosthetic clinics during the process of requesting advice or customer service or during the C-Leg prosthetists certification process.

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"Trial fittings are an appropriate and effective method to identify C-Leg responders. In the cohort of likely and preselected responders we observe a high degree of utilization of functional benefit. The degree of the utilization seems to be of similar magnitude when compared those observed in C-Leg populations. Responders can be found in all age groups, all mobility grades and all amputation etiologies.

Safety and specifically falls are a severe challenge in the transfemoral amputee population. The rates of subjects reporting annual falling is higher than reported before.

Age, mobility grade and amputation etiology do not correlate with a subjects capability to exhibit functional benefits when fitted with a C-Leg. Neither do any of these parameters possess relevant diagnostic power. We rather observe that trial fitting with a C-Leg does severely influences mobility grade rating itself. We therefore conclude that the assessment of an individual's potential is required to reveal the subjects potential and that none of those demographic factors and specifically the use of the evaluation of mobility grade rating are justified to deny such an assessment." (Hahn et al., 2015)

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© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.