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Highsmith MJ, Kahle JT, Miro RM, Cress ME, Lura DJ, Quillen WS.
School of Physical Therapy and Rehabilitation Sciences, University of
South Florida, Tampa, FL.

Functional performance differences between the Genium and C-Leg prosthetic knees and intact knees

Journal of Rehabilitation Research and Development 2016; 53(6): 753-766.

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
With Genium compared to C-Leg:

- **Genium scores significantly higher than C-Leg**
in the upper-body flexibility (UBF), balance (BAL) and endurance (END) domains
- **Genium scores higher than C-Leg**
upper-body and lower-body strength scores

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Genium scores higher than C-Leg in Physical Functional Performance

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UBS = upper-body strength; LBS = lower-body strength; UBF = upper-body flexibility; BAL = balance; END = endurance; ++ p < 0.05 (significant)

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
Interventional, randomized crossover design:

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“There were no significant differences in functional UBS between nonamputees and persons with TFA regardless of knee condition. Compared with the C-Leg, Genium use improved the UBF, BAL, and END domains of functional performance, likely because of improved confidence, willingness to lift and carry greater mass, and ability to move faster during activity. These benefits may be technologically due to the incorporation of a faster processing speed and axial load data assisting in regulating knee resistance and offering new functions such as stance locking and back-ward stepping. In the LBS, UBF, BAL, and END domains, C-Leg use resulted in significantly lower scores compared with nonamputees. Genium use significantly reduced the magnitude of impairment. The only domain in which persons with TFA performed significantly lower than nonamputees regardless of knee condition was the END domain. In terms of total CS-PFP10 performance, C-Leg use resulted in significantly lower function compared with nonamputees, whereas Genium use was not significantly different from nonamputees. Nonetheless, regardless of knee condition, persons with TFA did not equal or surpass nonamputees in any functional performance domain, suggesting room for improvements in TFA integrated functional performance. Further, the CS-PFP10 test was able to detect statistically significant differences of small effect size between prosthetic knee conditions, which should be interpreted with caution because the test has not been formally assessed in persons with TFA.” (Highsmith et al., 2016)

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