---

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Mileusnic MP, Hahn A, Reiter S.

Otto Bock Healthcare Products, Vienna, Austria.

## Effects of a Novel Microprocessor-Controlled Knee Kenevo on the Safety, Mobility, and Satisfaction of Lower-Activity Patients with Transfemoral Amputation

Journal of Prosthetics and Orthotics, 2017; vol. 29 (4):198-205

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**Kenevo**

---

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With Kenevo compared to previous fitting (mostly non-microprocessor controlled knees NMPKs)

→ **Reduced frequency of stumbles and falls**

Subjects that never stumble increased from 8% to 50%

Fear of falling reduced in 50% of subjects (negative trend)

→ **Wheelchair dependence reduced (from 87% to 37%)**

→ **Positive trend on Mobility and ADLS (Questionnaires, LCI-5, PLUS-M & Houghton scores)**

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
## Falls (N=12)

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
Stumbles and falls were evaluated by asking "how often do you stumble or fall with your current prosthesis?" at baseline measurements (with their old prosthesis) and again two months after being fit with Kenevo.

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Subjects:	29 (11 F) unilateral transfemoral amputees
Previous knee:	Polycentric (36%), brake (27%), locked (18%), MPK (18%) (N=22)
Amputation causes:	Vascular disease (46%), infection (38%), cancer (8%), trauma (8%) (N=13)
Mean age:	63.2 (± 9.5) years (N=16)
Mean time since amputation:	6.3 (± 8.9) years (N=12)
MFCL:	83% MFCL-2, 13% MFCL-3, 4% MFCL-1 (N=23)

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Demographic data and data on current prosthetic fitting was collected via questionnaire. Questions addressed topics such as pain, satisfaction and safety, together with validated clinical tests (Houghton, LCI-5, PLUS-M).

2 months after fitting, subjects returned for follow-up testing. The questions presented during baseline were again used to evaluate different aspects of Kenevo, together with questions comparing the subject's perception of Kenevo to their previous fitting.

Due to observation nature of this study, the completeness of data collected at different centres varied.

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"High patient satisfaction with Kenevo was observed in lower-activity individuals with transfemoral amputation. New knee functionalities designed specifically to target the safety challenges of this population seem to be effective and beneficial. The first clinical experiences collected during this observational study suggest that Kenevo offers several advantages to lower-activity persons with amputation, particularly in the areas of safety, functional mobility, preference, and satisfaction. These observations are in agreement with published evidence on effects of MPKs in limited community ambulators. Higher quality research is needed to confirm the conclusions reached in this observational study." (Mileusnic et al. 2017)

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© 2017, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.