## The Psychosocial and Biomechanical Assessment of Amputees Fitted with Commercial Multi-grip Prosthetic Hands – Case Study: Michelangelo hand

In Castelli, VP & Troncossi M, "Grasping the Future: Advances in Powered Upper Limb Prosthetics" (2012), Chapter 5, pp: 59-77. ISBN 978-1-60805-438-1

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With Michelangelo hand compared to Digital twin hand:


- **Michelangelo hand reduces compensatory movements**
- **Michelangelo hand gives more natural gesture and posture**
- **Patient is more satisfied with Michelangelo hand**

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### Disk task performance test with Michelangelo and Digital twin hand

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In disk task participant is moving the disk, positioned in front of the prosthetic hand, over the table. The participant moves disk over the table from the prosthetics hand to the sound hand, then in front of the participant and backwards. Michelangelo hand took same amount of time to perform disk task as sound hand (7s). Digital twin hand needed much more time to complete the same task (11s).

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Subjects: male, unilateral transradial amputee, dominant side  
 Previous: Digital Twin hand (Otto Bock)  
 Amputation causes: trauma  
 Mean age: 50 years  
 Mean time since amputation: 30 years

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"Results highlighted an increased satisfaction with the new multi-grip hand and, remarkably, the new prosthesis triggered a higher level of embodiment, with a mind changing in the use the previous hand as well. Thanks to pleasant appearance and functional features of Michelangelo, the patient started to assume more natural gestures and postures also with the traditional myoelectric hand, reporting this different way of thinking the prosthesis as "a fundamental step for an amputee". Regarding the biomechanical assessment, the shoulder biomechanics was positively influenced by the availability of the lateral grip and by the overall hand shape, which allowed the patient to approach cylindrical and coin-shaped objects in a more natural way, limiting the shoulder compensatory movements." (Cutti et al. 2012)

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